



SOUTHPAW®

Store these instructions in a safe place along with the enclosed maintenance checklist.

Instruction Sheet

Tower Ladder 1200 Ladder Wall 1250 Single Rung Climbing Ladder 1255



Take time to familiarize yourself with the use and maintenance of this equipment before using. Please file this Instruction Sheet in an accessible location: the Instruction Sheet must be read and familiar to anyone who is using Southpaw Enterprises Equipment.

IMPORTANT: Please check equipment immediately for any damage that may have resulted during shipping. If any damage has occurred, please contact Southpaw Enterprises, Inc. immediately.

Warning: Failure to properly install a Southpaw forged steel eyebolt (1/2" or 5/8") may result in serious injury or death. If you have any questions about the installation of the Southpaw forged steel eyebolt, please contact Southpaw Enterprises Customer Service at 1-800-228-1698.

MATERIALS LIST

- | | |
|---|---------------------------|
| 1 - Ladder | 1 - Instruction Sheet |
| 2 - Safety Snaps (1 - Snap w/Tower Ladder) | 1 - Maintenance Checklist |
| 2 - Mat straps w/Quick Links and foam covers attached | |

MINIMUM HEIGHT REQUIREMENT AND SPACE NEEDED

Minimum height required is 8 ft. Higher heights will require a height adjustment system. Assure a clear distance of 6 ft. from walls and other obstacles.

WORKING LOAD: 250 lbs.

SUSPENSION

Single and Tower Ladders require only a single point suspension system. Ladder Wall requires two suspension points.

- **Note:** Do not hang under light fixtures.
- **Special Note:** A Southpaw Safety Rotational Device (# 2000) should be used to prevent abnormal wear with single point suspension.

Assembly Procedures

- Attach the top safety snaps to ceiling eye bolt or suspension system.
- Position mat directly under ladder. Wrap adjustable mat straps around mat with attached quick links positioned directly under vertical ladder ropes. Tighten the straps snugly around mat. Attach the quick links to the bottom of the vertical ladder ropes. Slip foam covers from web down over quick links. These straps steady the ladder without making it rigid. When ladder is hanging properly, all rungs with the exception of the top and bottom rung will lift up. This allows the client to crawl through, to achieve a different climbing task. Before climbing occurs, mat straps should be placed around the width of the mat, placed far enough apart to maintain stability. Extra mats may be required to allow enough padding.
- **Special note:** Stay near to control movement of client. Watch for places that might pinch fingers.

MAINTENANCE CHECKLIST

All of the following wear points should be checked on a monthly basis and recorded on the checklist enclosed.

- Check ceiling eye bolt for wear and movement.
- Check safety snaps at top connection point for wear.
- Check mat straps.
- Check foam covers.
- Check ropes at top and bottom of bar.
- Check metal rope clips (the metal ring under each ladder bar).
- **Note:** If a rope clip comes loose or if a rope frays the ladder must be returned (for repair or replacement). If foam covers are damaged, please call Southpaw Enterprises and we will replace at no cost.

STORAGE PROCEDURE

- Store out of direct sunlight.
- Store where not accessible by clients without proper supervision.

CLEANING INSTRUCTIONS

- Dust or wipe with a damp cloth.
- Do not immerse in water.
- Dry before using.
- Do not wax ladder rungs.

- ◆ **Check ceiling eye bolt for wear and movement. Use a Southpaw Safety Rotational Device to prevent abnormal wear.**
- ◆ **Floor mats of proper thickness should be positioned properly under and around all suspended equipment.**

CAUTION

Safety Precautions - The safety of your clients is dependent on the appropriate choice of activity to match, yet challenge their physical abilities. Use appropriate spotting / safety techniques. In addition, the layout of your activity area, proper installation and maintenance of your equipment is critical. Southpaw equipment is not intended for recreational or outdoor use. Under no condition should this equipment be used without proper supervision.

This sheet is for example only. It is not meant to, in any way, replace local, county, state, or Federal standards. Be aware - it is your responsibility to find all pertinent safety information that pertains to your facility.

If you have any questions or comments, please call or write

Southpaw Enterprises, Inc.

P.O. Box 1047, Dayton, Ohio 45401-1047

Toll Free - (800) 228-1698

Please visit us at southpaw.com

© 2024 Southpaw Enterprises, Inc.