



Instruction Sheet

Tadpole Swing 121925

Store these instructions along with the enclosed maintenance checklist in a safe place. You may also access them on our website.



Take time to familiarize yourself with the use and maintenance of this equipment before using. Please file this Instruction Sheet in an accessible location. The Instruction Sheet must be read and familiar to anyone who is using Southpaw Enterprises Equipment.

IMPORTANT: Please check equipment immediately for any damage that may have resulted during shipping. If any damage has occurred, please contact Southpaw Enterprises, Inc. immediately.

Warning: Failure to properly install a forged steel eyebolt (1/2" or 5/8") may result in serious injury or death. If you have any questions about the installation of the forged steel eyebolt, please contact Southpaw Enterprises

MATERIALS LIST

- | | |
|------------------------------|-----------------------------------|
| 1 - Tadpole Swing | 1 - Instruction Sheet |
| 3 - 30 lb. Stimulation Loops | 1 - Safety Snap Instruction Sheet |
| 2 - Locking Safety Snaps | 1 - Maintenance Checklist |

MINIMUM HEIGHT REQUIREMENT AND SPACE NEEDED

- Minimum height requirement is 8 feet.
- Caution must be used to keep a clear distance of 6 feet from walls and other obstacles.

HEIGHT ADJUSTER AND SAFETY ROTATIONAL DEVICE NEEDED

- Use of the Southpaw Height Adjuster (#5000) is recommended for ease of adjusting your equipment.
- A Southpaw Safety Rotational Device (#2000) must be used for maximum effectiveness and decreased stress on equipment.

WORKING LOAD: 30 - 90 pounds (150 pounds with additional loops). You can add or remove stimulation loops to adjust the Tadpole Swing to the weight of the user to get the proper amount of vertical movement. Additional or replacement loops are available from Southpaw. Single loop, 30 lbs. # 1937 Double loop, 60 lbs # 1938

HANGING PROCEDURES

IMPORTANT: KEEP ALL PARTS TOGETHER

Attach your Tadpole Swing to your ceiling support system using the Locking Safety Snap provided. The Tadpole Swing is now ready for use.

MAINTENANCE CHECKLIST

All of the following wear points should be checked on a monthly basis and recorded on the checklist enclosed.

- Locking Safety Snaps - see enclosed Safety Snap Instruction Sheet
- Stimulation Loops - Check for loose stitching around the "ball" near the top and the bottom, and the stitches holding the web loops on the top and bottom. If loose or missing stitches are found, discontinue use until repaired or replaced.
- Check wood for cracks or splinters. If cracked, discontinue use immediately. For splinters, sand smooth and refinish with clear lacquer.
- Check all seams in fabric for loose or missing stitches. If loose or missing stitches are found, discontinue use until repaired.

MATTING PRECAUTIONS

- Crash mats and side mats must be used at all times.

CLEANING INSTRUCTIONS

- Check with your infectious disease department before using these procedures.
- Check for client allergies first.
- Remove the foam cushion and the nylon mesh fabric is machine washable / dryable (low heat only).

TO REMOVE THE SWING:

1. Open the hook and loop wrap that covers the connector strap (fig. 1).
2. Turn the hex nut to open the quick link fastener (fig. 2).
3. Remove the quick link fastener and slide the swing off of the strap (fig. 3).



Fig.1



Fig.2



Fig.3

TO REMOVE THE CUSHION:

1. Open the hook and loop fastener on the seat of the swing.
2. Reach inside and grasp the end of the cushion. Pull the end of the cushion out through the opening.

REASSEMBLY: Reassemble cushion and swing in reverse order.

STORAGE PROCEDURES

- Store out of direct sunlight.
- Store where not accessible without proper supervision.

CAUTION

- ◆ Check ceiling eye bolt for wear and movement. Use a Southpaw Safety Rotational Device to prevent abnormal wear.
- ◆ Floor mats of proper thickness should be positioned properly under and around all suspended equipment.

Safety Precautions - The safety of your clients is dependent on the appropriate choice of activity to match, yet challenge their physical abilities. Use appropriate spotting / safety techniques. In addition, the layout of your activity area, proper installation and maintenance of your equipment is critical. Southpaw equipment is not intended for recreational or outdoor use. Under no condition should this equipment be used without proper supervision.

This sheet is for example only. It is not meant to, in any way, replace local, county, state, or Federal standards. Be aware - it is your responsibility to find all pertinent safety information that pertains to your facility.

If you have any questions or comments, please call or write

Southpaw Enterprises, Inc.

P.O. Box 1047, Dayton, Ohio 45401 -1047

Toll Free - (800) 228-1698

Please visit us @ www.southpaw.com

© 2019 Southpaw Enterprises, Inc