



southpaw enterprises, inc.[®]

Store this instruction sheet in a safe place or you may access it on our website.

Instruction Sheet

Dolphin Air Pillow 1324

Walrus Air Mat 1325 - Seal Air Mat 1326



Take time to familiarize yourself with the use and maintenance of this equipment before using. Please file this Instruction Sheet in an accessible location: the Instruction Sheet must be read and familiar to anyone who is using Southpaw Enterprises Equipment.

IMPORTANT: Please check equipment immediately for any damage that may have resulted during shipping. If any damage has occurred, please contact Southpaw Enterprises, Inc. immediately.

MATERIALS LIST

1 - Air Mat
1 - Instruction Sheet

TOOLS NEEDED

Southpaw Inflator/Deflator or Vacuum Cleaner

MINIMUM SPACE NEEDED

Use care to assure a clear distance of six feet from walls and other obstacles. Under certain inflation, care must be used when bouncing so client does not hit low ceilings, light fixtures, etc. Use appropriate mats around air mats on floors and walls.

WORKING LOAD: 150 lbs.

Fig. A: Pull out cap and valve




Fig. B: Removal of Cap and Plug





Fig. C: Pushing cap and valve back into recess in Mat or Pillow



INFLATION PROCEDURES

Step 1. If cap and valve were pushed in during shipping, pull it out. (See figure A).

Step 2. Remove cap and plug. (See figure B). Inflate the Air Mat using Southpaw's Inflator/Deflator or a household vacuum cleaner. If using a household vacuum cleaner that has blower capabilities, make sure that the hose is clean (blowing it outside the room is advisable). This must be done before starting the operation. Inflate until desired firmness is obtained. Vary the amount of air in the mat to increase or decrease the mat's stability. Do not pressurize unit. Slightly over-inflate to allow for air loss while inserting the plug. Stop inflating at that point, drop the blower hose and put the plug in as quickly as possible. Hold with finger and start the cap on straight. Screw until tight. Then turn 1/8 turn more. If threads feel immediately tight, stop, unscrew cap completely, and restart.

NOTE: Do not over inflate. This is a non-elastic material. It will not maintain it's shape as it looses air like a rubber balloon. It is literally a bag of air. Over inflation will cause premature wear.

Step 3. When inflated, push internal cap back into the hole. (See figure C).

NOTE: Frequently Check seams and valves for obvious leaks. Repair or replace as necessary.

STORAGE PROCEDURES

- Store out of direct sunlight.
- Store away from all chemicals and household cleaners.
- Store where not accessible by clients without proper supervision.

CLEANING PROCEDURES

Check with your infectious disease department before using these procedures.

- Clean with a damp cloth using soap and water or disinfectant.
- Dry before using.

CAUTION

- ◆ Sharp objects should be removed from clothing and shoes before using mats or pillows.
- ◆ Do not fill with water.
- ◆ Do not use as a trampoline or crash mat.
- ◆ Do not jump on mattress from elevated structure such as ladder or table.
- ◆ Should be checked for proper inflation before every use.
- ◆ If used in water, mats and pillows are not adequate flotation.
- ◆ **MATTING PRECAUTIONS:** Always use appropriate matting. Air Mats and Pillows are bags of air and should be used on a surface with padding underneath them.
- ◆ If bag breaks or explodes in any way, client may come down abruptly and contact ground.
- ◆ Should not be used on concrete or hard floor unless padding is in place between floor and bag.

Safety Precautions - The safety of your clients is dependent on the appropriate choice of activity to match, yet challenge their physical abilities. Use appropriate spotting / safety techniques. In addition, the layout of your activity area, proper installation and maintenance of your equipment is critical. Southpaw equipment is not intended for recreational or outdoor use. Under no condition should this equipment be used without proper supervision.

This sheet is for example only. It is not meant to, in any way, replace local, county, state, or Federal standards. Be aware - it is your responsibility to find all pertinent safety information that pertains to your facility.



If you have any questions or comments, please call or write
Southpaw Enterprises, Inc.
P.O. Box 1047, Dayton, Ohio 45401 -1047
Toll Free - (800) 228-1698
Come see us on the World Wide Web at <http://southpawenterprises.com>