



southpaw enterprises, inc.®

Instruction Sheet Scooter Board XL 150110

Store these instructions along with the enclosed maintenance checklist in a safe place. You may also access them on our website.



Take time to familiarize yourself with the use and maintenance of this equipment before using. Please file this Instruction Sheet in an accessible location: the Instruction Sheet must be read and familiar to anyone who is using Southpaw Enterprises Equipment.

IMPORTANT: Please check equipment immediately for any damage that may have resulted during shipping. If any damage has occurred, please contact Southpaw Enterprises, Inc. immediately.

MATERIALS LIST:

1 - Scooter Board XL w/ Pull Strap	1 - Instruction Sheet
2 - Leg Wraps	1 - Maintenance Checklist

WORKING LOAD: 150 lbs.

NOTE: A pull loop is attached to the front underside of the Scooter Board so that the client may be pulled along by the pull strap provided.

LEG WRAPS:

On the bottom of the Scooter Board there are two strips of hook and loop fastener. The loop portion may be used as leg wraps when sitting on the Scooter Board or they can wrap around the body when laying down on the board. Attach by pressing one end of the loop strap onto the hook strip that is fastened to the bottom of the board and wrapping it over the rider's legs or body and hooking the other end of the strip on the opposite side of the board.

The Scooter Board XL will accommodate the Southpaw Scooter Board Backrest (# 1952) which also has positioning straps to secure the rider.

MAINTENANCE CHECKLIST:

All of the following wear points should be checked on a monthly basis and recorded on the checklist enclosed.

- Check tow strap for loose or missing stitches. If loose or missing stitches are found discontinue use.
- Check tow strap swivel hook to make sure it is functioning properly, if it is not discontinue use.
- Check cover for rips and exposed foam. If foam is showing, repair or discontinue use to prevent clients from ingesting foam.
- Check casters to make sure all are rolling properly. If casters are not rolling properly contact Southpaw.

MATTING PRECAUTIONS:

- Always use appropriate matting.

CAUTION:

- Eliminate obstructions that may be on the floor.
- No standing on the Scooter Board.
- Only one client at a time.
- Avoid riding stairs.
- Avoid riding backwards.

STORAGE PROCEDURES:

Store where not accessible by clients without proper supervision.

CLEANING INSTRUCTIONS:

- Check for client allergies first.
- Check with your infectious disease department before using these procedures.
- Clean with a damp cloth using soap and water or disinfectant. Dry before using.

Safety Precautions - The safety of your clients is dependent on the appropriate choice of activity to match, yet challenge their physical abilities. Use appropriate spotting / safety techniques. In addition, the layout of your activity area, proper installation and maintenance of your equipment is critical. Southpaw equipment is not intended for recreational or outdoor use. Under no condition should this equipment be used without proper supervision.

This sheet is for example only. It is not meant to, in any way, replace local, county, state, or Federal standards. Be aware - it is your responsibility to find all pertinent safety information that pertains to your facility.



Southpaw Enterprises Inc.®

If you have any questions or comments, please call or write
Southpaw Enterprises, Inc.
P.O. Box 1047, Dayton, Ohio 45401 -1047
Toll Free - (800) 228-1698
Come see us on the World Wide Web at <http://southpawenterprises.com>