



southpaw enterprises, inc.®

Store these instructions in a safe place or with the enclosed maintenance checklist

Instruction Sheet

Weighted Body Shaw - Aromatherapy 240050

Weighted Body Shaw - Unscented 240060



Take time to familiarize yourself with the use and maintenance of this equipment before using. Please file this Instruction Sheet in an accessible location: the Instruction Sheet must be read and familiar to anyone who is using Southpaw Enterprises Equipment.

IMPORTANT: Please check equipment immediately for any damage that may have resulted during shipping. If any damage has occurred, please contact Southpaw Enterprises, Inc. immediately.

MATERIALS LIST:

- 1 - Weighted Body Wrap
- 1 - Instruction Sheet
- 1 - Maintenance Checklist

CAUTION: Improper use of this product could result in burn or injury. Due to the difference in heating times of various microwaves, caution should be used when first heating Weighted Body Wrap until a safe amount of heat is achieved and the proper heating time for your microwave is determined. You may use towel layers between skin and Weighted Body Wrap if it is too warm. People suffering from nerve damage, circulation problems or paralysis should use the Weighted Body Wrap only as directed by a physician. The plastic storage bag for this product poses a suffocation hazard to small children.

INSTRUCTIONS: Fold shawl in half twice and heat in microwave on high for 2 minutes. May also be used as a cold compress. Place in plastic bag and then into freezer for 2-3 hours.

CLEANING INSTRUCTIONS:

Check with your infectious disease department before using these procedures.

- Check for client allergies first.
- Clean with a damp cloth using soap and water or disinfectant. Do not machine wash.
- Dry before using.

STORAGE PROCEDURES

- Store where not accessible by clients without proper supervision.

Safety Precautions - The safety of your clients is dependent on the appropriate choice of activity to match, yet challenge their physical abilities. Use appropriate spotting / safety techniques. In addition, the layout of your activity area, proper installation and maintenance of your equipment is critical. Southpaw equipment is not intended for recreational or outdoor use. Under no condition should this equipment be used without proper supervision.

This sheet is for example only. It is not meant to, in any way, replace local, county, state, or Federal standards. Be aware - it is your responsibility to find all pertinent safety information that pertains to your facility.



**Southpaw
Enterprises Inc.®**

If you have any questions or comments, please call or write
Southpaw Enterprises, Inc.

P.O. Box 1047, Dayton, Ohio 45401 -1047

Toll Free - (800) 228-1698

Come see us on the World Wide Web at <http://southpawenterprises.com>