



southpaw enterprises, inc.®

Instruction Sheet Fiber Optic Tunnel 303490

Store these instructions with the enclosed maintenance checklist in a safe place. You may also access them on our website.



Take time to familiarize yourself with the use and maintenance of this equipment before using. Please file this Instruction Sheet in an accessible location: the Instruction Sheet must be read and familiar to anyone who is using Southpaw Enterprises Equipment.

IMPORTANT: Please check equipment immediately for any damage that may have resulted during shipping. If any damage has occurred, please contact Southpaw Enterprises, Inc. immediately.

MATERIALS LIST:

1 - Fiber Optic Tunnel 1 - Instruction Sheet 1 - Maintenance Checklist

SPECIAL NOTES:

- Must be used with Fiber Optic Light Source - connect harness to light source.
- Use caution if using Tunnel with more than one client at a time.
- Do not allow clients to stand on barrel.
- Do not use with any ramp or platform.
- Although the Tunnel is made of foam, use caution that fingers do not get pinched where the barrel meets the floor.
- The foam used in the Tunnel is constructed of very sturdy foam, however if it is flattened or deformed severely, it may not recover its original shape.
- Fiber Optic Tunnel comes fully assembled

WORKING LOAD: 225 lbs.

MAINTENANCE CHECKLIST:

All of the following wear points should be checked on a monthly basis and recorded on the checklist enclosed.

- Check the vinyl cover for loose stitches, rips, or tears. If loose stitches, rips, or tears are found, discontinue use until repaired. CAUTION: Do not allow clients to pick foam from inside and ingest it.

STORAGE PROCEDURES:

- Store out of direct sunlight.
- It is best to store Fiber Optic Tunnel in an upright position.
- Do not store near chemicals or household cleaners.
- Store where not accessible by clients without proper supervision.

CLEANING INSTRUCTIONS:

- Check for client allergies first.
- Check with your infectious disease department before using these procedures.
- Wipe with a damp cloth using soap and water or disinfectant.

Safety Precautions - The safety of your clients is dependent on the appropriate choice of activity to match, yet challenge their physical abilities. Use appropriate spotting / safety techniques. In addition, the layout of your activity area, proper installation and maintenance of your equipment is critical. Southpaw equipment is not intended for recreational or outdoor use. Under no condition should this equipment be used without proper supervision.

This sheet is for example only. It is not meant to, in any way, replace local, county, state, or Federal standards. Be aware - it is your responsibility to find all pertinent safety information that pertains to your facility.



Southpaw Enterprises Inc.®

If you have any questions or comments, please call or write
Southpaw Enterprises, Inc.

P.O. Box 1047, Dayton, Ohio 45401 -1047

Toll Free - (800) 228-1698

Come see us on the World Wide Web at <http://southpawenterprises.com>