



Store these instructions with the enclosed maintenance checklist in a safe place. You may also access them on our website.

Instruction Sheet

Variable Axis Swing Footrest 4601



Take time to familiarize yourself with the use and maintenance of this equipment before using. Please file this Instruction Sheet in an accessible location: the Instruction Sheet must be read and familiar to anyone who is using Southpaw Enterprises Equipment.

IMPORTANT: Please check equipment immediately for any damage that may have resulted during shipping. If any damage has occurred, please contact Southpaw Enterprises, Inc. immediately.

MATERIALS LIST

- 1 - Variable Axis Swing Footrest
- 1 - Instruction Sheet
- 1 - Maintenance Checklist

ATTACHING THE FOOTREST

- With your Variable Axis Swing empty, tilt the seat back until it is parallel with the floor.
- Remove the rubber end caps from the long straight tubes of the footrest. These will be replaced later.
- Remove the footrest clamps from the end of the long straight tubes by loosening the knobs and separating the two halves. Do not confuse the clamp with the positioning plate. The positioning plates are round with many teeth.
- Place the clamps onto the inner hoop of the Variable Axis Swing (the hoop with the seat on it), just above the swing's positioning plates. The knobs should be facing away from the swing, and the clamp should extend downward. Do not tighten the clamps yet. At this point, the hole in the clamp that will accommodate the long straight footrest tubes should be below the seat. **Note: The footrest clamps may have numbers stamped into one end. The clamp halves with the same numbers must remain a pair, and the ends with the numbers must be kept together. Failure to do so may result in an improper fit of the clamp.**
- Place the footrest just under the seat of the Variable Axis Swing, with the zipper on the material covered end facing the floor.
- Slide the long straight tubes into the available hole in the clamp. Adjust the footrest to the desired length, and tighten the clamps. The clamp will be angled slightly inward.
- Replace the rubber end caps onto the ends of the long straight tubes of the footrest. At this point, the footrest should be firmly attached to the Variable Axis Swing.

CHANGING THE ANGLE OF THE FOOTREST

- Loosen the knobs on the positioning plates of the footrest until the teeth disengage. Change the angle of the cloth covered section to the desired position, and tighten the knobs until the teeth fully engage. **Do not loosen the knobs more than necessary to adjust the footrest.**

CHANGING THE LENGTH OF THE FOOTREST

- To adjust how close the material covered end of the footrest is to the seat, loosen the knobs on the footrest clamps (not the positioning plates) and slide the long straight tubes of the footrest to the desired position. Tighten the clamps until the footrest is once again firmly held by the clamps. **Loosen the clamps only enough to slide the tubes to the proper position.**

MAINTENANCE CHECKLIST

All of the following wear points should be checked on a monthly basis and recorded on the checklist enclosed. Check cover for torn or ripped seams. If tears or rips are found, discontinue use until repaired.

STORAGE PROCEDURES

- Store out of direct sunlight.
- Store where not accessible by client without proper supervision.

CLEANING INSTRUCTIONS

Check with your infectious disease department before using these procedures.

- Check for client allergies first.
- Cover is machine washable - cold water, mild detergent.
- Tumble dry, low heat.

Safety Precautions - The safety of your clients is dependent on the appropriate choice of activity to match, yet challenge their physical abilities. Use appropriate spotting / safety techniques. In addition, the layout of your activity area, proper installation and maintenance of your equipment is critical. Southpaw equipment is not intended for recreational or outdoor use. Under no condition should this equipment be used without proper supervision.

This sheet is for example only. It is not meant to, in any way, replace local, county, state, or Federal standards. Be aware - it is your responsibility to find all pertinent safety information that pertains to your facility.

If you have any questions or comments, please call or write

Southpaw Enterprises, Inc.

P.O. Box 1047, Dayton, Ohio 45401-1047

Toll Free - (800) 228-1698

Please visit us at southpaw.com

© 2024 Southpaw Enterprises, Inc.