



# Instruction Sheet

## ASwing Frame 7290

Store these instructions with the enclosed maintenance checklist in a safe place. You may also access them on our website.



Take time to familiarize yourself with the use and maintenance of this equipment before using. Please file this Instruction Sheet in an accessible location: the Instruction Sheet must be read and familiar to anyone who is using Southpaw Enterprises Equipment.

**IMPORTANT: Please check equipment immediately for any damage that may have resulted during shipping. If any damage has occurred, please contact Southpaw Enterprises, Inc. immediately.**

The goal of the ITN Frame is to offer the therapist the ability to use suspended therapy in places not possible because of lack of a permanent suspension point.

Another goal is to give the itinerant therapist quick and easy set up and take down, thus allowing the therapist more time with the client and less time assembling equipment.

### MATERIALS LIST

- 1 - ITN Frame
- 1 - Instruction Sheet
- 1 - Maintenance Checklist

### SETTING UP THE FRAME

- Spread all four legs as far as allowed by the cables at the bottom of the legs. Hang equipment from the eye bolt.

### CAUTION

- Do not swing beyond the angle of the legs of the frame. Failure to comply may cause the frame to tip over, resulting in serious injury to the client and / or the therapist.
- Do not hang from or hang anything on the legs of the frame. It is not designed to support weight in this manner.
- Be sure to conduct swinging activity so that the path of the swing is between the frame legs.
- Do not allow the client to push off from or grab onto the legs while swinging. This could lead to loss of control and possible injury.

## TAKING DOWN THE FRAME

- Swing the leg in toward the center of the frame and lay the frame down. Collapse the three remaining legs.

## TO REMOVE THE TOP COVER FOR INSPECTION

- Unsnap the four snaps located on the underside of the top plate. Lift the cover off of the frame.

## MAINTENANCE CHECKLIST

All of the following wear points should be checked on a weekly basis and recorded on the checklist enclosed.

- ✓ Remove the top cover and inspect the swivel eye bolt. Check that the nut holding the eye bolt has not loosened, and the eye bolt turns freely. If the nut has loosened or the eye bolt does not turn freely, discontinue use immediately and contact Southpaw for repairs.
- ✓ Check eye bolt for wear. If wear is more than 30% of the diameter of the eye bolt, discontinue use and contact Southpaw for replacement.
- ✓ Check the eight copper stop swages to be sure that the cable is secure. If the cable has slipped in the swage, discontinue use immediately and contact Southpaw for repairs.
- ✓ Check that the cables at the bottom of the legs are secure. If cable is not secure, discontinue use immediately until repaired.

**WORKING LOAD:** 300 lbs.

**CAUTION:** When transporting the ITN frame, you should take care to set the frame down *on its side*. If you need to set it on one of its ends, please do so **GENTLY**. Dropping the frame down on either end repeatedly will cause damage to the legs, or to the cap. By taking the time to set the frame down on its side, you will greatly extend the life of the frame and reduce the chances of needing early repairs.

## CAUTION

- ◆ **Floor mats of proper thickness should be positioned properly under and around all suspended equipment.**

**Safety Precautions** - The safety of your clients is dependent on the appropriate choice of activity to match, yet challenge their physical abilities. Use appropriate spotting / safety techniques. In addition, the layout of your activity area, proper installation and maintenance of your equipment is critical. Southpaw equipment is not intended for recreational or outdoor use. Under no condition should this equipment be used without proper supervision.

**This sheet is for example only.** It is not meant to, in any way, replace local, county, state, or Federal standards. Be aware - it is your responsibility to find all pertinent safety information that pertains to your facility.

If you have any questions or comments, please call or write  
Southpaw Enterprises, Inc.  
P.O. Box 1047, Dayton, Ohio 45401 -1047  
Toll Free - (800) 228-1698  
Please visit us @ [www.southpaw.com](http://www.southpaw.com)  
© 2019 Southpaw Enterprises, Inc.