



# southpaw enterprises, inc.®

Store these instructions along with the enclosed maintenance checklist in a safe place. You may also access them on our website.

## Instruction Sheet

### *AdvantageLine* Scooter Board 150050



Take time to familiarize yourself with the use and maintenance of this equipment before using. Please file this Instruction Sheet in an accessible location: the Instruction Sheet must be read and familiar to anyone who is using Southpaw Enterprises Equipment. If this Instruction Sheet is misplaced, please call or write for a replacement. (See address below).

**IMPORTANT:** Please check equipment immediately for any damage that may have resulted during shipping. If any damage has occurred, please contact Southpaw Enterprises, Inc. immediately.

#### MATERIALS LIST:

- 1 - Scooter Board w/ Pull Strap
- 1 - Instruction Sheet
- 1 - Maintenance Checklist

**WORKING LOAD:** 150 lbs.

**NOTE:** A pull loop is attached to the front underside of the Scooter Board so that the client may be pulled along by the pull strap provided.

#### MAINTENANCE CHECKLIST:

All of the following wear points should be checked on a monthly basis and recorded on the checklist enclosed.

- Check tow strap for loose or missing stitches. If loose or missing stitches are found discontinue use.
- Check tow strap swivel hook to make sure it is functioning properly, if it is not discontinue use.
- Check cover for rips and exposed foam. If foam is showing, repair or discontinue use to prevent clients from ingesting foam.
- Check casters to make sure all are rolling properly. If casters are not rolling properly contact Southpaw.

#### MATTING PRECAUTIONS:

- Always use appropriate matting.

#### STORAGE PROCEDURES:

- Store where not accessible by clients without proper supervision.

#### CLEANING INSTRUCTIONS:

- Check for client allergies first.
- Check with your infectious disease department before using these procedures.
- Clean with a damp cloth using soap and water or disinfectant. Dry before using.

#### CAUTION:

- Eliminate obstructions that may be on the floor.
- No standing on the Scooter Board.
- Only one client at a time.
- Avoid riding stairs.
- Avoid riding backwards.

**Safety Precautions** - The safety of your clients is dependent on the appropriate choice of activity to match, yet challenge their physical abilities. Use appropriate spotting / safety techniques. In addition, the layout of your activity area, proper installation and maintenance of your equipment is critical. Southpaw equipment is not intended for recreational or outdoor use. Under no condition should this equipment be used without proper supervision.

**This sheet is for example only.** It is not meant to, in any way, replace local, county, state, or Federal standards. Be aware - it is your responsibility to find all pertinent safety information that pertains to your facility.



If you have any questions or comments, please call or write  
Southpaw Enterprises, Inc.  
P.O. Box 1047, Dayton, Ohio 45401 -1047  
Toll Free - (800) 228-1698  
Come see us on the World Wide Web at <http://southpawenterprises.com>