



southpaw enterprises, inc.[®]

Store these instructions along with the enclosed maintenance checklist in a safe place. You may also access them on our website.

Instruction Sheet: Bean Chair 156420



Take time to familiarize yourself with the use and maintenance of this equipment before using. Please file this Instruction Sheet in an accessible location. The Instruction Sheet must be read and familiar to anyone who is using Southpaw Enterprises Equipment.

IMPORTANT: Please check equipment immediately for any damage that may have resulted during shipping. If any damage has occurred, please contact Southpaw Enterprises, Inc. immediately.

MATERIALS LIST:

- | | |
|------------------------------|---------------------------|
| 1 - Bean Chair (unassembled) | 1 - Fabric Seat |
| 2 - Side Boards | 1 - Allen Key |
| 3 - Dowels | 1 - Instruction Sheet |
| 6 - Screws | 1 - Maintenance Checklist |

ASSEMBLY:

Select one side of the Bean Chair. Place one end of a dowel into the pocket on the side board. From the opposite side, insert a screw into the hole and using the allen key provided, tighten the screw. Place remaining dowels into the next two pockets on the same side and secure with a screw. While the side board is in rocking position, slide the longer piece of fabric over the top two dowel rods. The remaining fabric will slide over remaining dowel. Take the other side of the bean chair, align dowels into pockets and secure with screws.

WORKING LOAD: 200 lbs

INSTRUCTIONS FOR USE:

Chair may be used on either side. One side will provide a stable base, the other will provide a rocking motion.

MAINTENANCE CHECKLIST:

All of the following wear points should be checked on a monthly basis and recorded on the maintenance checklist enclosed.

- Check wood for cracks or splinters. If cracked, discontinue use. For splinters, sand until smooth and refinish with clear lacquer.
- Check for loose screws. These may need to be tightened periodically, but do not over tighten.
- Check for torn cover or ripped seams. If found, discontinue use until repaired.

MATTING PRECAUTIONS:

Always use appropriate matting.

STORAGE PROCEDURES:

- Store where not accessible by clients without proper supervision.

CLEANING INSTRUCTIONS:

Check with your infectious disease department before using these procedures.

- Check for client allergies first.
- Wipe wood with a damp cloth using soap and water or disinfectant. Dry before using.
- Cover is machine washable - Cold water, mild detergent. Air dry - no heat.

Safety Precautions - The safety of your clients is dependent on the appropriate choice of activity to match, yet challenge their physical abilities. Use appropriate spotting / safety techniques. In addition, the layout of your activity area, proper installation and maintenance of your equipment is critical. Southpaw equipment is not intended for recreational or outdoor use. Under no condition should this equipment be used without proper supervision.

This sheet is for example only. It is not meant to, in any way, replace local, county, state, or Federal standards. Be aware - it is your responsibility to find all pertinent safety information that pertains to your facility.



**Southpaw
Enterprises Inc.®**

If you have any questions or comments, please call or write
Southpaw Enterprises, Inc.

P.O. Box 1047, Dayton, Ohio 45401 -1047

Toll Free - (800) 228-1698

Come see us on the World Wide Web at <http://southpawenterprises.com>