



southpaw enterprises, inc.®

Store these instructions in a safe place or with the enclosed maintenance checklist

Instruction Sheet

Fine Motor Activity Table 191000



Take time to familiarize yourself with the use and maintenance of this equipment before using. Please file this Instruction Sheet in an accessible location: the Instruction Sheet must be read and familiar to anyone who is using Southpaw Enterprises Equipment.

IMPORTANT: Please check equipment immediately for any damage that may have resulted during shipping. If any damage has occurred, please contact Southpaw Enterprises, Inc. immediately.

MATERIALS LIST:

1 - Activity Table
8 - Storage Containers w/ Lids

1 - Instruction Sheet
1 - Maintenance Checklist

NO ASSEMBLY REQUIRED

- Table top height is adjustable from 18" to 26" in 2" increments. Adjust the height, one leg at a time, by loosening the two knobs holding each leg. Slide the knobs into the vertical slot then into the horizontal slot at the height of your choice. Tighten both knobs. Adjust all four legs in the same manner.
- The table has holders for storage bins on one side. Lids must be removed to fit into storage racks.
- The table has a 1" lip around the edges to prevent objects from rolling off table top.
- The table has a cut out on one side to provide access to table top.
- The table has a cut out on one end and is predrilled for optional ramp (#151000).

MAINTENANCE CHECKLIST:

All of the following wear points should be checked on a monthly basis and recorded on the checklist provided.

- Check wood for cracks or splinters. If cracks are found, discontinue use until repaired. Sand splinters smooth and refinish with clear lacquer.
- Check that all knobs and screw are tight. If knobs or screws are found loose, discontinue use until tightened.

CLEANING INSTRUCTIONS:

Check with your infectious disease department before using these procedures.

- Check for client allergies first.
- Wipe with damp cloth using soap and water or disinfectant. Dry before using.

Safety Precautions - The safety of your clients is dependent on the appropriate choice of activity to match, yet challenge their physical abilities. Use appropriate spotting / safety techniques. In addition, the layout of your activity area, proper installation and maintenance of your equipment is critical. Southpaw equipment is not intended for recreational or outdoor use. Under no condition should this equipment be used without proper supervision.

This sheet is for example only. It is not meant to, in any way, replace local, county, state, or Federal standards. Be aware - it is your responsibility to find all pertinent safety information that pertains to your facility.



**Southpaw
Enterprises Inc.®**

If you have any questions or comments, please call or write
Southpaw Enterprises, Inc.
P.O. Box 1047, Dayton, Ohio 45401 -1047
Toll Free - (800) 228-1698
Come see us on the World Wide Web at <http://southpawenterprises.com>



southpaw enterprises, inc.®

Store these instructions in a safe place or with the enclosed maintenance checklist

Instruction Sheet

Ramp - Fine Motor Activity Table 151000



Take time to familiarize yourself with the use and maintenance of this equipment before using. Please file this Instruction Sheet in an accessible location: the Instruction Sheet must be read and familiar to anyone who is using Southpaw Enterprises Equipment.

IMPORTANT: Please check equipment immediately for any damage that may have resulted during shipping. If any damage has occurred, please contact Southpaw Enterprises, Inc. immediately.

MATERIALS LIST:

- | | |
|-----------------------------------|---------------------------|
| 1 - Ramp | 1 - Instruction Sheet |
| 2 - 5/16-18 x 1/2" Hex Head Bolts | 1 - Maintenance Checklist |
| 2 - 5/16" Split Lock Washers | |

TO ATTACH THE RAMP:

- Locate the two attachment holes on the end of the ramp (the end with the cut out). Each hole has a threaded insert to screw the bolt into.
- Position the ramp smooth side up, with the hinge end of the ramp next to the mounting holes.
- Put one lock washer on each of the hex head bolts.
- Raise the hinge end of the ramp. With the hinge turned down, put a bolt through the hole in one end of the hinge and start the bolt into the threaded insert in the table. Repeat on the opposite side of the ramp.
- Tighten both bolts securely using a 1/2" or an adjustable wrench.

MAINTENANCE CHECKLIST:

All of the following wear points should be checked on a monthly basis and recorded on the checklist provided.

- Check wood for cracks or splinters. If cracks are found, discontinue use until repaired. Sand splinters smooth and refinish with clear lacquer.
- Check that all knobs and screw are tight. If knobs or screws are found loose, discontinue use until tightened.

CLEANING INSTRUCTIONS:

Check with your infectious disease department before using these procedures.

- Check for client allergies first.
- Wipe with damp cloth using soap and water or disinfectant. Dry before using.

Safety Precautions - The safety of your clients is dependent on the appropriate choice of activity to match, yet challenge their physical abilities. Use appropriate spotting / safety techniques. In addition, the layout of your activity area, proper installation and maintenance of your equipment is critical. Southpaw equipment is not intended for recreational or outdoor use. Under no condition should this equipment be used without proper supervision.

This sheet is for example only. It is not meant to, in any way, replace local, county, state, or Federal standards. Be aware - it is your responsibility to find all pertinent safety information that pertains to your facility.



**Southpaw
Enterprises Inc.®**

If you have any questions or comments, please call or write
Southpaw Enterprises, Inc.

P.O. Box 1047, Dayton, Ohio 45401 -1047

Toll Free - (800) 228-1698

Come see us on the World Wide Web at <http://southpawenterprises.com>