

Store these instructions along with the enclosed maintenance checklist in a safe place. You may also access them on our website.

Instruction Sheet In-FUN-ity Ladder Climbing Wall 152020



Take time to familiarize yourself with the use and maintenance of this equipment before using. Please file this Instruction Sheet in an accessible location: the Instruction Sheet must be read and familiar to anyone who is using Southpaw Enterprises Equipment.

IMPORTANT: Please check equipment immediately for any damage that may have resulted during shipping. If any damage has occurred, please contact Southpaw Enterprises, Inc. immediately.

MATERIALS LIST:

1 - Vertical Ladder (2 pieces)

1 - 7/32" Hex Wrench-XH90150

2 - Connector Boards

1 - Instruction Sheet

8 - Button Head Cap Screws 3/8-16 x 1-1/2"-XH90420 1 - Maintenance Checklist

8 - 3/8" Flat Washers-XH7034

ASSEMBLY:

- Put a 3/8" flat washer on four of the 3/8" x 1-1/2" cap screws. Position the connector board over the joint between the top and bottom side rails with the straight edge toward the back with the holes in the board aligned with the screw inserts in the side rails as shown in fig.1. Put one of the 3/8" bolts into each of the holes in the connector board and screw them into the corresponding screw insert in the side rails. Tighten securely using the hex key wrench provided.
- Assemble the rails and connector board to the opposite side in the same manner.

MAINTENANCE CHECKLIST:

All of the following wear points should be checked on a monthly basis and recorded on the checklist enclosed.

- Check that all screws and bolts are tight. If any screws or bolts are found to be loose, discontinue use until they are tightened.
- Check wood surfaces for cracks and splinters. If cracks are found, discontinue use immediately and contact Southpaw Enterprises. Splinters may be sanded smooth and refinished with clear lacquer.

WORKING LOAD: 300 lbs

MATTING PRECAUTIONS:

Always use appropriate matting under and around equipment.

STORAGE PROCEDURE:

Make sure that equipment is not accessible by clients without proper supervision.

CLEANING INSTRUCTIONS:

- Check with your infectious disease department before using these procedures.
- Check for client allergies first.
- Clean with a damp cloth using soap and water or disinfectant.
- Dry before using.

Figure 1



Safety Precautions - The safety of your clients is dependent on the appropriate choice of activity to match, yet challenge their physical abilities. Use appropriate spotting / safety techniques. In addition, the layout of your activity area, proper installation and maintenance of your equipment is critical. Southpaw equipment is not intended for recreational or outdoor use. Under no condition should this equipment be used without proper supervision.

This sheet is for example only. It is not meant to, in any way, replace local, county, state, or Federal standards. Be aware - it is your responsibility to find all pertinent safety information that pertains to your facility.