



Store these instructions in a safe place along with the enclosed maintenance checklist. You may also access them on our website.

# Instruction Sheet

## Padded Frog Swing Seat 1927



Take time to familiarize yourself with the use and maintenance of this equipment before using. Please file this Instruction Sheet in an accessible location: the Instruction Sheet must be read and familiar to anyone who is using Southpaw Enterprises Equipment.

**IMPORTANT:** Please check equipment immediately for any damage that may have resulted during shipping. If any damage has occurred, please contact Southpaw Enterprises, Inc. immediately.

### MATERIALS LIST

- 1 - Padded Frog Swing Seat
- 1 - Instruction Sheet
- 1 - Maintenance Checklist

### INSTRUCTIONS FOR USE

Open the hook and loop fastener and place the padded seat centered on top of the frog swing seat. Fold the flap with the hook fastener under the seat, then fold the flap with the loop (fuzzy) fastener onto the hook.

### MAINTENANCE CHECKLIST

All of the following wear points should be checked on a monthly basis and recorded on the checklist enclosed.

- Check for loose stitches or ripped seams. If loose stitches or ripped seams are found, discontinue use until repaired.

### STORAGE PROCEDURES

- Store out of direct sunlight.
- Store where not accessible by clients without proper supervision.

### CLEANING INSTRUCTIONS

- Check with your infectious disease department before using these procedures.
- Check for client allergies first.
- Wipe with a damp cloth using soap and water of disinfectant.
- Dry before using.

**Safety Precautions** - The safety of your clients is dependent on the appropriate choice of activity to match, yet challenge their physical abilities. Use appropriate spotting / safety techniques. In addition, the layout of your activity area, proper installation and maintenance of your equipment is critical. Southpaw equipment is not intended for recreational or outdoor use. Under no condition should this equipment be used without proper supervision. **This sheet is for example only.** It is not meant to, in any way, replace local, county, state, or Federal standards. Be aware - it is your responsibility to find all pertinent safety information that pertains to your facility.

If you have any questions or comments, please call or write  
Southpaw Enterprises, Inc.  
P.O. Box 1047, Dayton, Ohio 45401-1047 • Toll Free - (800) 228-1698  
Please visit us at [southpaw.com](http://southpaw.com)  
© 2024 Southpaw Enterprises, Inc.

