



# southpaw enterprises, inc.®

Store these instructions along with the enclosed maintenance checklist in a safe place. You may also access them on our website.

## Instruction Sheet

## Vibro-Acoustic Podium #303112



Take time to familiarize yourself with the use and maintenance of this equipment before using. Please file this Instruction Sheet in an accessible location. The Instruction Sheet must be read and familiar to anyone who is using Southpaw Enterprises Equipment.

**IMPORTANT: Please check equipment immediately for any damage that may have resulted during shipping. If any damage has occurred, please contact Southpaw Enterprises, Inc. immediately.**

### **MATERIALS LIST**

- Vibro-Acoustic Podium
- Podium Pad
- 24' speaker wire

### **INSTRUCTIONS FOR ASSEMBLY**

The Podium comes out of the box completely assembled. In order to make the Podium vibrate the product must be connected to a Southpaw Stereo Rover, Southpaw Corner or Wall Radio . Insert the RCA plug into jack on the Southpaw Stereo Rover or Radio. Use with any other stereo equipment will void the warranty.

### **CLEANING / INSPECTION / MAINTENANCE**

Recommended cleaning solution is a 50/50 mixture of Hydrogen Peroxide and water and/or a small amount of mild detergent. Use a soft cotton cloth in conjunction with any cleaner.


### **MAINTENANCE**

There is no regular maintenance necessary other than keeping the product clean. If vinyl should tear or rip, contact Southpaw to replace cover.

**CARE SHOULD BE TAKEN TO REMOVE ALL SHARP OBJECTS FROM POCKETS, AS WELL AS ANY METAL BRACES, PRIOR TO SITTING ON THE PODIUM.**

**Safety Precautions** - The safety of your clients is dependent on the appropriate choice of activity to match, yet challenge their physical abilities. Use appropriate spotting / safety techniques. In addition, the layout of your activity area, proper installation and maintenance of your equipment is critical. Southpaw equipment is not intended for recreational or outdoor use. Under no condition should this equipment be used without proper supervision.

**This sheet is for example only.** It is not meant to, in any way, replace local, county, state, or Federal standards. Be aware - it is your responsibility to find all pertinent safety information that pertains to your facility.



**Southpaw  
Enterprises Inc.®**

If you have any questions or comments, please call or write  
Southpaw Enterprises, Inc.

P.O. Box 1047, Dayton, Ohio 45401 -1047  
Toll Free - (800) 228-1698

Come see us on the World Wide Web at  
<http://southpawenterprises.com>